

# SHARESOURCE CONNECTIVITY PLATFORM for use with **HOMECHOICE CLARIA** APD SYSTEM

## NEW **SHARESOURCE** CLINICAL PORTAL ENHANCEMENTS



### IMPROVED ULTRAFILTRATION (UF) CLARITY & REPORT EFFICIENCIES

**Sharesource**, by Baxter, is designed to offer a higher standard of care for patients on home dialysis.

By applying innovative technology, **Sharesource** creates a connected home PD therapy experience to make a meaningful difference in our patients' lives.

The latest version of **Sharesource** brings enhanced functionality with greater UF data clarity and report generation time-savings.



#### Ultrafiltration (UF) Insight

A complete look at UF. Maintaining adequate fluid balance is critical to improve patient outcomes<sup>1</sup>.



#### Reporting Efficiency

Data comprehension made easy. The right data in the right spot. With 1-click, create a ready to use organized and neatly formatted report.

1. European Best Practice Guidelines; NDT Vol 2, suppl 9, Dec 2005, 1-37.

# Sharesource offering a higher standard of care for patients on home dialysis

Patients using APD devices with **Sharesource** had trends indicating fewer hospitalizations than patients using APD devices without **Sharesource**<sup>1</sup>.

1. Firanek et. al. 2017; Disclaimer: This is an observational study of 189 patients using APD devices with **Sharesource**.



## Ultrafiltration (UF) Complete Picture

Get a comprehensive look at your patient's UF including initial drain UF and manual exchange UF. Methodically packaged to easily interpret results. This UF detail is also included in the corresponding Treatment Summary Report.

Day Therapy UF	
Manual Exchange UF	450 ml
Initial Drain UF	1400 ml
<b>Total Day UF</b>	<b>1850 ml</b>
Night Therapy UF	
<b>Total Night Cycle UF</b>	<b>1150 ml</b>
Total Reported Therapy UF	
<b>Total UF</b>	<b>3000 ml</b>

UF reported is based onycler recorded and manually entered data. The Clinician should ensure any unreported IIF data is accounted for



## Prescribed vs. Actual Smarter Reporting

As a key tool to understanding your patient's treatment, the Prescribed vs. Actual report has been enhanced to include night UF, total UF and pre-weight/blood pressure vitals. Also, reformatted to a shorter length, averaged per period and revised column titles for easier interpretation.

**Wellness Clinic: Homechoice Claria Prescribed Versus Actual Report**

Date	Device Program Name	Night Therapy Time (H:MM)	Total Therapy Volume (ml)	Avg Fill Volume (ml)	No. of Night Cycles	Avg Dwell Time (H:MM)	Initial Drain Volume (ml)	Night UF/Total UF (ml)	Pre-Weight (kg)	Pre-Sys/Dia (mm/Hg)
22 Apr 2018	Rx: Weekend	8:30	14,000	2000	6	1:40	2200	--/--	Target: 86	--
22 Apr 2018	Weekend	8:25	14,050	1950	6	1:35	2505	1750/2200	88.2	138/83
23 Apr 2018	Rx: Daily	7:00	12,000	2000	5	1:14	2100	--/--	Target: 86	--
23 Apr 2018	Daily	7:05	12,075	1940	5	1:10	1950	1690/1700	87.7	139/84
23 Apr 2018	Daily	6:55	12,050	1930	5	1:02	2000	1980/2050	88.1	140/85
23 Apr 2018	Daily	7:08	12,100	1960	5	1:05	2050	1611/1700	87.7	139/84
<b>Period Average</b>	--	<b>7:23</b>	<b>12,569</b>	<b>1945</b>	<b>5.25</b>	<b>1:13</b>	<b>2126</b>	<b>1758/1913</b>	<b>87.9</b>	<b>139/84</b>



## Effortless Reporting

Burden-free and efficient reporting with 1-click of a button. The Clinical and Patient Treatment History reports are organized and pre-formatted saving you time.

1	Clinic:	Global Sharesource	Report Tir	04 May 2018, 12:05
2	Report Na	Homechoice Claria P	Exported	Garcia, Christy
3	Patient N:	Branham, Kristi	Patient Confidential Information	
4	Attending Dr.	Ken, Doctor	Report Stc	#####
5	Date of Bl	9-Jul-88	Report En	2-May-18
6	Clinic Patient ID:			
7	Baxter Pat O-	TRAIN-1002071		
9	Date	Device Pr	Night The	Total Ther
10	#####	Weekday	5:02	10,002
11	#####	Weekday	4:41	10,007



1	Clinic:	The OTTAWIA Hospital_CaseOf	Report Timestamp	04 May 2018, 17:05
2	Report Name:	Homechoice Claria Patient Treatment Hi	Exported By	Garcia, Christy
3	Patient Name:	Admas, John	Patient Confidential Information	
4	Attending Healthcar	JONES, MITTON	Report Start Date	17 April 2018
5	Date of Birth:	02 March 1960	Report End Date	02 May 2018
6	Clinic Patient ID:			
7	Baxter Patient ID:	O-ST123-1002618		
10	Date	Device Programme Name	Night Therapy Time (H:MM)	Total Therapy Volume (ml)
11	20 April 2018	DP1	0:00	0
12	21 April 2018	DP1	1:00	798
13	21 April 2018	DPME	0:48	799

Clinicians switching from APD to APD with **Sharesource** spent a greater proportion of time on proactive patient care<sup>2</sup>.

2. Bunch et. al. 2017; Disclaimer: Some bias may exist due to study design; no statistical difference was demonstrated for technical failure.

