Medium cut-off dialyzer improves erythropoiesis stimulating agent resistance in a hepcidin-independent manner in maintenance hemodialysis patients : results from a randomized controlled trial

Lim JH, Jeon Y, Yook JM, et al. Sci Rep. 2020 Sep 29;10(1):16062.

Enabled by

Theranova

Medium cut-off dialyzer improves erythropoiesis stimulating agent resistance in a hepcidin-independent manner in maintenance hemodialysis patients : results from a randomized controlled trial

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Post-hoc analysis of the prospective, randomized, controlled, open-label trial to compare the ESA resistance between MCO and high-flux HD patients. The original randomized controlled trial was designed to evaluate the effectiveness of MCO dialyzer on quality of life compared to high-flux dialyzer. Chronic maintenance HD patients treated with high-flux dialyzer at the Kyungpook National University Hospital were enrolled from July 2018 and followed up for 12 weeks. The MCO group changed their dialysis membrane from a high-flux dialyzer to a MCO dialyzer (Theranova 400). The control group continued their treatments with the high-flux dialyzer.

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

Comparison of the ERI difference¹

Serial changes in the ERI¹

The p values for the difference between groups from baseline to 12 weeks were calculated in the generalized estimating equation models. Asterisk (*) indicates significant difference between groups at the time point ($p < 0.05$).

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

Reduction ratio of serum hepcidin and TNF- α ¹

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The iron metabolism regulatory pathway¹

Red arrows indicate dominant effects. *It improves the iron status by both hepcidin-independent and dependent pathway, although independent pathway was more dominant in our study.

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ERI=erythropoietin resistance; ESA=erythropoiesis stimulating agent; ESRD=end-stage renal disease; HD=hemodialysis; HF=high-flux; MCO=medium cut off; OL-HDF=online hemodiafiltration; RR=reduction ratio; TNF-α=tumor necrosis factor-alpha. TSAT=transferrin saturation
 Reference. 1. Lim JH, Jeon Y, Yook JM, Choi SY, Jung HY, Choi JY, Park SH, Kim CD, Kim YL, Cho JH. Medium cut-off dialyzer improves erythropoiesis stimulating agent resistance in a hepcidin-independent manner in maintenance hemodialysis patients: results from a randomized controlled trial. Sci Rep. 2020 Sep 29;10(1):16062.