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Integrated Treatment System
: HDx-HC treatment option Vablet

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ated Treatment System
HC treatment option Vablet

» ARTIS PHYSIO SYSTEM

Individualized treatment with the Artis Physio™ dialysis system

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 Brochure
Download

» ARTIS PHYSIO SYSTEM

Individualized treatment with the Artis Physio™ dialysis system

Artis Physio™ 특징

Individualized Quality-assured Dialysis (IQAD)

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Individualized treatment: overcoming the hemodialysis challenges

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Making possible personal.

ARTISET

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Baxter
Novaline
BLOODLINE TUBING SYSTEM

BLOODLINE
TUBING SYSTEM
SIMPLE
EFFICIENT

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


Hemodialysis

**HEMODIALYSIS DOUBLE NEEDLE
ARTISET HD DNL HC**

The ARTISET HD DNL HC blood tubing system for double needle hemodialysis treatment includes:

- The HEMOSCAN cassette to measure blood volume variations, enabling the **HemoCentral** modality
- A connector on the venous bypass line for ULTRA Prime function
- A pre-connected infusion line for saline priming

ARTISET blood tubing system for HD treatment	PRODUCT DESCRIPTION	ULTRA HD LINE
Product code	86278	86283
Volume (ml)	100	100
Shelf life	2 years	2 years
Dimensions	See data sheet	See data sheet
Units per box	25	25

Materials

Body	PC (polypropylene) medical grade, USP Class II	Body	PC (polypropylene) medical grade, USP Class II
Connectors	PC (polypropylene) medical grade, USP Class II	Connectors	PC (polypropylene) medical grade, USP Class II
Injection site	Latex free rubber	Injection site	Latex free rubber

On-line HDF – HF in Post dilution

**HEMODIALYSIS DOUBLE NEEDLE
HDF POSTDILUTION
ARTISET HD DNL HC**

The ARTISET HD DNL HC can also be used for HDF postdilution including the **UltraControl** modality.

Accessory

ON-LINE SUBSTITUTION FLUID ULTRA HDF LINE

The ULTRA HDF LINE includes a single-use ultrafilter which removes bacteria from fluid and ensures a sterile path from its membrane to the point of infusion. Furthermore, it provides fresh sites for adsorption of endotoxins.

This bloodline accessory is dedicated for on-line metabolites, to be used in combination with the ARTISET range and includes:

- A single-use ultrafilter for production of substitution fluid
- A cassette body for auto-loading operation
- The T connector for rinse back to facilitate the rinse back operation

ARTISET blood tubing system for HDF and HF treatment	PRODUCT DESCRIPTION	ULTRA HDF LINE
Product code	86278	86283
Volume (ml)	100	100
Shelf life	2 years	2 years
Dimensions	See data sheet	See data sheet
Units per box	25	25

Materials

Body	PC (polypropylene) medical grade, USP Class II	Body	PC (polypropylene) medical grade, USP Class II
Connectors	PC (polypropylene) medical grade, USP Class II	Connectors	PC (polypropylene) medical grade, USP Class II
Injection site	Latex free rubber	Injection site	Latex free rubber

Perforances

Max. Water (ml)	100L
-----------------	------

 **Brochure**

 **Download**

BLOODLINE TUBING SYSTEM

ARTISET

NOVALINE


On-line HDF – HF in Pre/Post dilution

**HEMODIALYSIS DOUBLE NEEDLE
HDF PRE AND POSTDILUTION
HF PRE-DILUTION
ARTISET HDF PREPOST**

The ARTISET PREPOST blood tubing system for double-needle hemodialysis, HDF pre and postdilution, or HF pre-dilution treatment includes:

- A connector on the arterial dialyzer line for on-line pre-dilution
- A connector on the venous dialyzer line for ULTRA Prime function and HDF postdilution, including the **UltraControl** modality
- The **HEMOCONTROL** cassette to measure blood volume variations, enabling the **HemoControl** modality

Accessory

ON-LINE SUBSTITUTION FLUID ULTRA HDF LINE

The ULTRA HDF LINE includes a single-use ultrafilter which removes bacteria from fluid and ensures a sterile path from its reservoir to the point of infusion. Furthermore, it provides fresh sites for adsorption of endotoxins.

This bloodline accessory is dedicated for on-line modalities, to be used in combination with the ARTISET range and includes:

- A single-use ultrafilter for production of substitution fluid
- A cassette body for auto-bleeding operation
- The T-connector for rinse back to facilitate the rinse-back operation

ARTISET blood tubing system for HDF and HF treatment		Accessory	
Product code	Shelf life	Product code	Shelf life
85277	3 years	15213	1 year
85278	3 years	15214	1 year
85279	3 years	15215	1 year
85280	3 years	15216	1 year

Materials		Materials	
Tubing	PVC polyurethane (medical grade), 30MPa tensile	Filter	Membrane: PES/PTMDS Polypropylene
Cassette body	PVC polyurethane (medical grade), 30MPa tensile	Frame	Polycarbonate
Injection site	Latex free rubber	Sealing	PVC polyurethane (medical grade), 30MPa tensile

Performances	
Max. infusion rate	100 L/h

INTEGRATED BLOOD CIRCUIT FOR SIMPLICITY – IN JUST ONE MOVE

As a pioneer and leading innovator in dialysis therapy, we present the ARTISET blood tubing system that brings new ways to raise treatment quality, staff satisfaction and operational efficiency compared to conventional blood tubing systems.

This ARTISET blood tubing system is an exclusive blood circuit designed to work with the **Artis Physio** system. Compared with conventional bloodlines, the compact and highly integrated design offers various important advantages both from a medical and handling perspective.

ARTISET blood tubing system

HD DN – ARTISET blood tubing system		Accessory	
Product code	Shelf life	Product code	Shelf life
85278	3 years	15213	1 year
85279	3 years	15214	1 year
85280	3 years	15215	1 year
85281	3 years	15216	1 year

On-line HDF/HF in Pre/Post dilution – ARTISET blood tubing system		Accessory	
Product code	Shelf life	Product code	Shelf life
85277	3 years	15213	1 year
85278	3 years	15214	1 year
85279	3 years	15215	1 year
85280	3 years	15216	1 year

On-line HDF in Post dilution – ARTISET blood tubing system		Accessory	
Product code	Shelf life	Product code	Shelf life
85278	3 years	15213	1 year
85279	3 years	15214	1 year
85280	3 years	15215	1 year
85281	3 years	15216	1 year

Note: A leading group check all safety requirements recommended in case of TDS (Free backflow control test).

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

Download

» BLOODLINE TUBING SYSTEM

Baxter
Novaline
BLOODLINE TUBING SYSTEM

BLOODLINE TUBING SYSTEM

SIMPLE

EFFICIENT

NOVALINE MODELS FOR HD ON AK SYSTEMS

CODE	MODEL	SYSTEM	DRIVEN BY BOX	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT LINE	INFUSION SET + CONNECTOR	TRANSFUSION PRE-FILTER	PRE-DILUTION CHALLENGE	POST-DILUTION	VENOUS DILUTION	DIALYZER CONNECTOR	REMOVABLE LABELING LINE
Y50170	SL 050	AK 91 S, AK 94, AK 95	DI	24	Etoam	GDHP Free	3x12 mm	No	No	1	Yes	1	Dual	No
Y50172	SL 91 DI	AK 91 S, AK 94, AK 95	DI	24	Etoam	GDHP Free	3x12 mm	No	No	1	Yes	1	Dual	No

NOVALINE MODELS FOR HD, HF AND HDF ON AK SYSTEMS

CODE	MODEL	SYSTEM	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT LINE	INFUSION SET + CONNECTOR	TRANSFUSION PRE-FILTER	PRE-DILUTION CHALLENGE	POST-DILUTION	VENOUS DILUTION	DIALYZER CONNECTOR	REMOVABLE LABELING LINE	
Y50173	SL 91 DI V1	AK 900 S, AK 900 Ultra S	24	Etoam	GDHP Free	3x12 mm	No	No	1	Yes	Yes	2	Dual	No

For safe and proper use of products mentioned herein, please refer to the appropriate Instructions for Use.

Accessories: Novaline
CE 0086 **CE 0197**

Baxter, AK 91 S, AK 94, AK 95, AK 900 S, AK 900 Ultra S, AK 91, Aklix, Physis, Integra and Novaline are trademarks of Baxter International Inc. or its subsidiaries.

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SPECIALIZED FOR
저체중[LOW BODY WEIGHT] 환자

Polyflux 2H/6H

Dialyzer Portfolio Brochure Download

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1. Schepers E, Doreau G, Essel S, et al. Assessment of the association between increasing membrane pore size and endotoxin permeability using a novel experimental dialysis simulation set up. BMC Nephrology. 2016; 19:1.
 2. Baxter. Evaluation of synthetic membranes for blood purification: the case of the Polyflux Family. Regional Dial Transplant 1993; 18(Suppl 7):10-18.
 3. Baxter. Data on file. Skyline Sales Report. 2016.
 4. Collin-Bernard CE, et al. Cytotoxic effects exerted by polysulfone dialyzer membranes depend on different sterilization processes. Int J Artif Organs. 2015; 43:483-495.
 5. O'Anderson PR, et al. Polyamide dialyzer allergy in dialysis patients. Regener Dial. 1977; 2:144-148.
 6. Kizawa S, et al. Polyamide Membrane for Medical Applications. Chemtech International. 2003; 10:125-132.
 7. Talemans C, et al. Clinical assessment of Performance and Blood Compatibility Profile of a New Synthetic Low Flux Hemodialyzer. Blood Purif. 2010; 29:214-215.

The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (PFI) quality assurance system of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/45/EC.
 For safe and proper use of the device, please refer to the Instructions for Use
CE 0086

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Baxter
Polyflux H

Polyflux H

DESIGNED FOR: **CONNECTIVE (HDF-HF)** | OTHER APPLICABLE THERAPIES: **HFHD (High Flux HD)** | MEMBRANE: **POLYAMIX (PAES/PVP/PA, BPA-free)**

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* According to EN 1083/ISO 8537
 * UF-Coefficient: measured with bovine blood, wet 32%, PCl at 37°C
 * β₂-m: measured at Q_d 150 mL/min, Q_d 150 mL/min, UF at 100 mL/min
 * Setting coefficients: measured with bovine for human β₂m, Q_d 300 mL/min, UF at 100 mL/min
 * Creatinine in vitro: measured at Q_d 150 mL/min, wet 32%
 * HDI/HTI: made: measured at UF=40 mL/min, Q_d 150 mL/min
 1. Kwon C, et al. Evaluation of synthetic membranes for blood purification: the case of the Polyflux Hemofiltration Dialyzer. *Transfusion* 2003;43(10):1492-99.
 2. Malhotra P, et al. High efficiency proteinase enzyme hemofiltration reduces all-cause mortality in hemodialysis patients. *J Am Soc Nephrol*. 2012;23:448-491.
 3. Fisher J et al. Are all dialyzers compatible with the convective volumes suggested for modification in the hemodialysis literature? *Am J Kidney Dis* 2014;64:95-104.
 4. Baxter. *Baxter on file*. *Baxter's Safety Report*. 2010.
 5. Calkins-Bonnie EE, et al. Optimal efficacy of polypropylene dialyzer membranes depend on different sterilization processes. *Int J Artif Organs*. 2011; 34:485-490.
 6. Karamantzi PP, et al. Endotoxin reduction in dialysis patients. *Blood Purif* 1993;13:161-163.
 7. Schepers E, Oltmanns G, Elert S, et al. Assessment of the association between increasing membrane pore size and protein permeability using a novel experimental dialyzer simulation set up. *BMC Nephrology*. 2016; 17:1.
 8. Gruchis V, et al. Effect of CO₂ clearance on uremic toxins and clearance on beta₂-m in low hemofiltration treatments. *Int J Artif Organs* 2012; 35:343-345.
 9. Meert N, et al. Effective removal of protein-bound paraprotein by different convective modalities: comparison of high-flux dialyzer. *Transfusion* 2005; 45:442-449.

The products meet the applicable provisions of Annex I (Essential Requirement) and Annex II (Quality assurance system) of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/45/EC.
 For safe and proper use of the device, please refer to the Instructions for Use.
CE 0086

Baxter

서울시 용문구 용문1로 100 10018 Tel: 02. 6262. 7100 Fax: 02. 6272. 7101 www.baxter.co.kr

Polyflux
Capillary Dialyzer
210H

04/2017/000

» BAXTER DIALYZER PORTFOLIO

THE POLYFLUX DIALYZER

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Baxter
Polyflux 2H/6H
DIALYZER

Polyflux 2H/6H

DESIGNED FOR: HFHD [High flux HD] | OTHER APPLICABLE THERAPIES: CONVECTIVE [HDF-HF] | MEMBRANE: POLYAMIX [PAES/PVP, PA, BPA-free]

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1. Sano C, et al. Evaluation of synthetic membranes for blood purification in the case of the Polyflux family. *Regener Dial Transplant* 2003; 35:105-108.

2. Goldstein GS, et al. Polyflux® dialyzer: a new option for small children requiring dialysis for ESRD. *Peritoneal Dialysis* 2003; 35:371-376.

3. Chertow GM, et al. Pediatric Dialysis, Chapter 19. Maintenance hemodialysis during infancy. *2012 Kidney International Supplements*, Chapter 19, 25-46.

4. Kessler DM, et al. A new synthetic dialyzer with advanced permeability for enhanced low-molecular-weight protein removal. *Am J Nephrol* 2008; 30:541-554.

5. Liu L, et al. Barrier function of low and high flux synthetic membranes for endotoxins. *Contaminated Dialysis Unit: Blood Purif* 2003; 25:206.

6. Christensen FF, et al. Ethylene oxide sterilization in dialysis patients. *Regener Dial* 1973; 3:141-143.

The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (EU quality assurance system of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/45/EC)

For safe and proper use of the device, please refer to the Instructions for Use
CE 0086

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* According to EN 13813:2015 B.2.1.
 † UF-Coefficient: measured with bovine blood, 300 mL, Prt 300, 20°C
 ‡ K_{0.45} urea: calculated at Q_B=300 mL/min, Q_D=500 mL/min, UF=40 mL/min
 § Seiving coefficients: measured with bovine plasma, Q_B=300 mL/min, UF=40 mL/min
 ¶ Clearances in vitro: measured at UF=40 mL/min, 45°C

1. Baxter, REVACLEAR White Paper, USMFM02140002, May 2015.
 2. Moore A, et al. Poster SPW1, presented at 15th ERA-EDTA congress, Istanbul (Turkey), 2015.
 3. Kuo C, et al. Evaluation of synthetic membrane for blood purification: the case of the Prisma family. Regional Dial Transplant 2013;18(5):71-78.
 4. Ward S, et al. Abstract SA-0010, presented at the 43rd ASN congress, San Francisco (USA), 2017.
 5. Sheng B, et al. Effect of increasing dialyzer flow rate on ultrafiltration volume, phosphate and beta2-microglobulin during clinical hemodialysis. Regional Dial Transplant 2015; 20:990-995.
 6. Yu G, et al. Information as a cause of malnutrition, albuminemia cardiovascular disease, and poor outcome in hemodialysis patients. Hemodial Int 2014; 8:108-121.
 7. Baxter. Data on file. Biocompatible water use calculation, 2015.
 8. Baxter. REVACLEAR dialyzer priming guide. 30470112_C_2016

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vidences on HDx

oL outcomes
JH, et al. Sci Rep. 2020)

HDx : ENABLED BY THE THERANOVA DIALYZER

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* According to EN 12815/ISO 8037
 UF-Coefficient: measured with bovine blood, at 32% PC, pH 7.3, 37°C
 K_v area: calculated at Q_b 300 mL/min, Q_d 100 mL/min, UF-CO 50 mL/min
 Sieving coefficients: measured with human plasma, Q_b 300 mL/min, UF-CO 50 mL/min
 Clearances in vitro: measured at UF-CO 50 mL/min, 40% to 20% CrCl₂, C_{in} 0.8% NaCl

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1. Baines C, et al. The role of Expanded Hemodialysis. Blood Purif 2015; 41:1-11
 2. Hocher CM, et al. The Rationale for Expanded Hemodialysis. Therap Apher Dial. 2015; 19:143-52
 3. Kitzler A, et al. Performance of hemodialysis with novel medium cut-off dialyzer. Nephrol Dial Transpl 2015; 30:1165-73
 4. Buechler M, Kitzler A, et al. MCO membranes: enhanced laboratory analysis by Case Sciences. Reports 2015; 5:1844
 5. Zenger L, et al. Medium cut-off membranes - closer to the natural kidney removal function. Int J Artif Organs 2015; 40(7):229-34
 6. Baxter. Data on file. Theranova Limited Confidential Disclosure Report 2014
 7. Baxter. Theranova 400/500 Instructions For Use. NDD-148-Rev 003, 2017-10-29
 8. Schepers E, Gijzen G, Elvel S, et al. Assessment of the association between increasing membrane pore size and reduction permeability using a novel experimental dialysis stimulation set-up. BMC Nephrology 2016; 17:3
 9. Cantalupo G, et al. Removal of large-molecule antibodies on expanded hemodialysis (EHD): a multicenter observational study of a months follow-up. ASN 2018 Kidney Week Abstract TH-PO027
 10. Kishimasa M, et al. Trial evaluating mid cut-off value membrane clearance of albumin and light chains in hemodialysis patients (MEMOVAL-FDS): a safety and efficacy study. ASN 2018 Kidney Week Abstract TH-PO023
 11. Savelkoul M, et al. Quality of life reported by patients with expanded hemodialysis by the Theranova dialyzer at RIZ Columbia. ASN 2018 Kidney Week Abstract TH-PO216

The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (P_u) quality assurance system of the Council Directive 93/42/EEC of 14 June 1993, amended by Directive 2007/47/EC

2020-01-29-2019

» HDx : ENABLED BY THE THERANOVA DIALYZER

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Ref. 1. Hocher A, et al. Performance of hemodialysis with novel medium cut-off dialysis. *Nephrol Dial Transplant* 2017; 32(1):165-72. 2. Cantalupo V, et al. Survival of large middle molecules on expanded hemodialysis. *Kidney Int* 2016; 89(5):1011-1018. 3. Cantalupo V, et al. Quality of hemodialysis and patient-reported outcomes on expanded hemodialysis. *Theranova dialyzer* (2015). Cambridge, MA: Baxter Healthcare Corporation. 4. Molyneux R, et al. The efficacy of off-line urea reduction ratio of dialysis on patient-reported outcomes. *SEMPER* (2015). A study of off-line study. *Am J Kidney Dis* 2015; 65(5):783-790. 5. Nesterov A, et al. The cost-effectiveness of hemodialysis versus hemodiafiltration in the Connective Tissue Study. *Nephrol Dial Transplant* 2016; 31(10):1872-1879.

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Baxter
U9000
ULTRAFILTER

U9000 ULTRAFILTER

Designed to produce ultrapure dialysis fluid*

* Intended to be used in conjunction with a water treatment system

DIACLEAR

U9000 ULTRAFILTER

INDICATIONS

The U9000 ultrafilter is indicated for purification of incoming water for dialysis fluid as well as purification of dialysis fluid to obtain microbiological high quality dialysis fluid. The risk of exposure to bacteria and endotoxins can thereby be minimized.

The U9000 ultrafilter is intended to be used in conjunction with a water treatment system.

MATERIALS

Membrane	PAES/PVP
Filtering material	Polyurethane (PUR)
Housing and caps	Polycarbonate (PC)
Gasket	Silicone rubber

CHARACTERISTICS

Effective membrane area (m ²)	2.4
Priming volume lumen (ml)	1.5
Priming volume filtrate side (ml)	210
Compatibility	The U9000 ultrafilter is only compatible with the Baxter HD systems equipped with the U9000 holder.

TYPICAL RETENTION VALUES FOR BACTERIA AND ENDOTOXINS

Challenge	LRV*
Bacterial challenge: <i>Pseudomonas dimorpha</i> ATCC 17148	> 7
Endotoxin challenge: E. coli 055:B5 endotoxin (Wako, USA)	> 3.5

* LRV = log10 number of organisms in challenge suspension/number of organisms in filtrate

PERFORMANCES

Typical pressure drop

Q_v = 500 ml/min - dialysis fluid at temperature = 27°C

~ 100 mmHg

Lifetime

On the AK10 V2, V3 and Artis dialysis systems:
The U9000 ultrafilter shall be changed after 3 months of usage or 150 disinfection cycles, whichever comes first.
On the AS 14 and AS 16/17 dialysis systems:
The U9000 ultrafilter shall be changed after 2 months of usage or 100 disinfection cycles, whichever comes first.
Maximum 12 disinfections with sodium hypochlorite over the lifetime of the U9000.
Maximum 8 disinfections with sodium carbonate over the lifetime of the U9000.

Disinfection / deactivation

The U9000 ultrafilter is disinfected with the dialysis machine. Follow the disinfection program recommended in the operator manual of the dialysis machine.

The U9000 ultrafilter is compatible with the following chemical agents:

- Peracetic acid (0.3, 1%)
- Sodium carbonate (0.5%)
- Sodium hypochlorite (0.5%)
- Citric acid (0.25%)

The U9000 ultrafilter is compatible with heat disinfection.

PRODUCT INFORMATION

Product code	Designation	Unit per box
11252	U9000	16
11604	U9000 AFAC	16

For safe and proper use of the device, please refer to the Instructions for Use

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Baxter, AK 14, AS 16, Artis and U9000 are trademarks of Baxter International, Inc. or its subsidiaries.

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Baxter
Diaclear
ULTRAFILTER

DIACLEAR ULTRAFILTER

Designed to produce ultrapure dialysis fluid*

* Intended to be used in conjunction with a water treatment system

DIACLEAR

DIACLEAR ULTRAFILTER

INDICATIONS

The **Diaclear** ultrafilter is indicated for purification of incoming water for dialysis fluid as well as purification of dialysis fluid to obtain microbiological high quality dialysis fluid. The risk of exposure to bacteria and endotoxins can thereby be minimized.

The **Diaclear** ultrafilter is intended to be used in conjunction with a water treatment system.

MATERIALS

Membrane	PSES
Plating material	Polyurethane (PU)
Housing	Polycarbonate (PC)

CHARACTERISTICS

Effective membrane area [m ²]	1.3
Priming volume (sum) [ml]	98
Priming volume filtrate side [ml]	135
Sterilization	Gamma rays

Compatibility: The **Diaclear** ultrafilter is only compatible with the Baxter HD system equipped with the **Diaclear** holder.

TYPICAL RETENTION VALUES FOR BACTERIA AND ENDOTOXINS

	LRV*
Bacterial challenge	> 7
Pseudomonas diminuta ATCC 19146	> 7
Endotoxin challenge: E. coli O157:H5 endotoxin (Whittaker, USA)	> 3.5

* LRV = \log_{10} number of organisms in challenge suspension / number of organisms in filtrate

PERFORMANCE

Typical pressure drop

Q _v = 500 ml/min - dialysis fluid at temperature = 27°C	< 100 mmHg
--	------------

Lifetime

The **Diaclear** ultrafilter shall be changed routinely after 30 dialysis sessions or after maximum 100 heat disinfections, depending on the results from the microbiological controls of the ultrapure dialysis fluid.

Disinfection / decalcification

The **Diaclear** ultrafilter is disinfected with the dialysis machine. Follow the disinfection program recommended in the operator manual of the dialysis machine.

The **Diaclear** ultrafilter is compatible with the following chemical agents:

- Peracetic acid (0.1%)
- Acetic acid (0.1%)
- Sodium hypochlorite (0.5%)
- Citric acid (0.2%)

The **Diaclear** ultrafilter is compatible with heat disinfection.

PRODUCT INFORMATION

Product code	Designation	Unit per box
104987	Diaclear Ultrafilter	12

For safe and proper use of the device, please refer to the Instructions for Use

CE 2797

Baxter and Diaclear are trademarks of Baxter International Inc. or its subsidiaries.

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Gambro
cleancart

CARTRIDGES

With **CleanCart A** and **CleanCart C** cartridges, regular dialysis machine maintenance and disinfection assure proper dialysis machine function and requires less intensive machine operation from clinical staff.

- Less space for storage and less heavy lifting for staff compared to liquid disinfectant
- Minimizes risk of personal injury and material damage when handling chemical agents
- Easy to streamline disinfection combining heat disinfection with decalcification and cleaning
- Removes organic deposits, fats and proteins with weekly use of **CleanCart A** cartridge¹
- Removes precipitated calcium and magnesium salts with daily use of **CleanCart C** cartridge¹
- Recyclable polypropylene container for less environmental impact compared to the liquid disinfectant

Dialysis machine maintenance and disinfection becomes a convenient and simple operation, so you can make the most of your time with patients.

Baxter

1. Operator manual

Baxter, Gambro and CleanCart are trademarks of Baxter International Inc. or its subsidiaries. 1000-00270004

HEMOCLEAN

Gambro **CleanCart** cartridge – specifications

CLEANCART A
The **CleanCart A** cartridge contains 13 grams of anhydrous sodium carbonate powder. The **CleanCart A** cartridge must be used in combination with the heat disinfection program to clean the dialysis machine fluid path. When circulated throughout the fluid path, the sodium carbonate solution removes organic deposits, fats and proteins¹.

CLEANCART C
The **CleanCart C** cartridge contains 32 grams of citric acid anhydrous powder. The **CleanCart C** cartridge must be used in combination with the heat disinfection program to both decalcify and disinfect the dialysis machine fluid path simultaneously. When circulated throughout the fluid path, the citric acid solution removes precipitated salts of calcium and magnesium¹.

TECHNICAL DATA		
	CleanCart A	CleanCart C
Action in combination with the heat disinfection program	Cleaning of the fluid path from organic deposits, fats, and proteins	Disinfection & decalcification of the fluid path
Cartridge content	13 grams	32 grams
Total cartridge weight	Approx. 32 grams	Approx. 50 grams
Cartridge and label material	Recyclable polypropylene (PP)	Recyclable polypropylene (PP)
Recommended usage	Refer to operator's manual for corresponding dialysis machine	Refer to operator's manual for corresponding dialysis machine
pH (prepared solution)	Approx. 11	Approx. 2
Residual test needed	No	No
Biodegradable content	Yes	Yes
Storage conditions	Store below +30°C	Store below +30°C
Shelf life	24 months from date of manufacture	24 months from date of manufacture
Product code	114012	114013
Packaging	Box of 33 cartridges	Box of 33 cartridges
Box dimensions	279 x 247 x 81 mm	279 x 247 x 81 mm
Total box weight	1200 grams	1800 grams

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Baxter, Gambro and CleanCart are trademarks of Baxter International Inc. or its subsidiaries. 1000-00270004

» WATER SYSTEM

BACK TO OVERVIEW

WRO 300H

Baxter

WRO300H

WATER PURITY FOR SINGLE PATIENTS

Automated procedures minimize patient and caregiver exposures
 Dialysis machines can start and stop the Gambro WRO 300H in situ
 Automatically shut-off at completion of disinfection process
 The WRO 300H unit can be programmed to automatically start heat disinfection at specified intervals


- Hot water disinfection
- Low noise level

Making possible personal.

[Brochure](#)

[Download](#)

CWP 64, 66

Baxter

Gambro

CWP60

DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}


- User-friendly control panel
- Simple, reliable and compact design
- Integration into Gambro dialysis machines
- Low noise level
- The CWP 60 system is available in five capacities

Making possible personal.

[Brochure](#)

[Download](#)

CWP 106H

Baxter

Gambro

CWP100

DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

The Gambro CWP 100 system is an integral part of the Hygienic Chain concept and allows for achieving and maintaining of the fluid purity.³

- Automated heat disinfection of the distribution loop
- Fully enclosed central water plant
- Water saving functionality
- Chemical disinfection of the RO membranes
- Quiet, reliable operation

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
WRO300H

WATER PURITY FOR SINGLE PATIENTS

- Automated procedures minimize patient and caregiver involvement
- Dialysis machine can start and stop the Gambro WRO 300 H unit
- Automatic shut-off at completion of disinfection process
- The WRO 300 H unit can be programmed to automatically start heat disinfection at specified intervals
- Hot water disinfection
- Low noise level*

Making possible personal.

The WRO 300 H unit

WATER FOR DIALYSIS
The quality of the water used in the preparation of dialysis fluid is very important. Even water considered as acceptable according to existing tap water regulations may have chronic as well as acute effects on the dialysis patient.^{1,2} The Gambro single patient reverse osmosis monitor WRO 300 H is designed to provide the high quality water needed for dialysis.^{2,3}

REVERSE OSMOSIS
Reverse osmosis is today the preferred method for the purification of water for dialysis. This method removes more than 96% of dissolved salts and more than 99% of all particles, bacteria and pyrogens in the water. Most tap waters can therefore be purified to a standard, which complies with existing recommendations for water for dialysis.^{2,3}

TECHNICAL DATA

PRODUCT WATER	
Output	Maximum 1.1 liter at +10°C and 0.15 MPa (1.5 bar) water pressure
Quality	Depends on inlet water quality. If potable water is used, and WRO 300 H is maintained according to the manual, the following minimum rejection rates will be obtained: Total dissolved salts: > 96% Bacteria and pyrogens: > 99%

FEED WATER SUPPLY	
Input	Min. 3.0 liter
Pressure	150-800 kPa (1.5-8 bar)
Temperature	+ 5 to + 20°C
Quality	Potable water should be used. Softener followed by carbon/particle filter ensures optimum performance. To insure maximum membrane life expectancy the following levels should not be exceeded:
Hardness	< 0.37 dH @ ppm as CaCl ₂
Iron	< 0.1 mg/l
Manganese	< 0.1 mg/l
Jackson Turbidity Unit (JTU)	< 1.2 TU
Total dissolved solids (TDS)	< 1000 mg/l
Silica density index (SDI)	< 5
Chlorine (free)	< 0.1 mg/l

DRAIN REQUIREMENTS	
Operation	1.2 g/l 1 liter
Peak flow (times)	Min. 3.0 liter required

CONNECTIONS	
Supply and drain lines	Designed for flexible, reinforced tubing, 8.0 mm x 2.3 mm
Product water loop	Designed for flexible, reinforced tubing, 8.0 mm x 3 mm

Specifications subject to change without prior notice.

For further information and operating instructions, please refer to applicable operator's manual.

1. Operator manual
2. Henschel et al., The Importance of water quality and contamination field connection, Blood Purification, 2006, 24, 11-18
3. ISO 13629:2014

CE 0086 This product is CE marked in accordance with the requirements in EC Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

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The WRO 300 H unit is a reverse osmosis unit designed specifically for dialysis. It combines simplicity, reliability and ease of use and is based on the long time experience of water treatment equipment within Gambro.

INTEGRATED HEAT OR CHEMICAL DISINFECTION
When a WRO 300 H unit is fitted to a Gambro dialysis machine, one of the disinfection programs will allow an integrated heat, alternatively chemical disinfection of the reverse osmosis unit, the connection line to the dialysis machine and the dialysis machine itself. This "end-to-end" action will ensure that the hygienic chain remains unbroken.

USER INTERFACE DISPLAY	
Product water conductivity	Temperature compensated product water conductivity, operating range 1-500 µS/cm
Feed water conductivity	Temperature compensated feed water conductivity, operating range 15-2000 µS/cm
Rejection rate	Rejection rate, operating range 0-100%
Flow	Flow and time, total run time, flow and heat disinfection, cleaning, etc.

TEMPERATURE MEASUREMENT	
Operating range	0-135°C

REVERSE OSMOSIS MEMBRANE	
Material	Polyamide, thin film composite
Configuration	Spiral-wound
g/salience	2:1:1

DISINFECTION & CLEANING	
Heat disinfection	User initiated or automatic start of heat disinfection
Chemical disinfection	Automatic initiation of disinfection. Rinse memory forcing the rinse program to start after chemical disinfection
Cleaning	Customised programs for different models

POWER SUPPLY	
Rated voltage	100, 115 or 220-240 V ± 10%, 50 or 60 Hz
Power	220-240 V: max 1800 W 100 V: max 1500 W 115 V: max 1500 W

AMBIENT	
Temperature	+ 5 to + 40 °C

DIMENSIONS	
Depth	Min: 520 mm Footprint: 380 mm
Width	Min: 240 mm Footprint: 495 mm
Height	565 mm
Weight	30 kg

0000-00000000

» WATER SYSTEM

WRO 300H

Baxter

WRO300H
WATER PURITY FOR SINGLE PATIENTS

Automated procedures minimize patient and caregiver exposures
Dialysis machines can start and stop the Gambro WRO 300H in situ
Automatic shut-off at completion of disinfection process
The WRO 300H unit can be programmed to automatically start heat disinfection at specified intervals

Hot water disinfection
Low noise level

Making possible personal.

[Brochure](#)

[Download](#)

CWP 64, 66

Baxter

Gambro
CWP60
DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

- User-friendly control panel
- Simple, reliable and compact design
- Integration into Gambro dialysis machines
- Low noise level
- The CWP 60 system is available in five capacities*

Making possible personal.

[Brochure](#)

[Download](#)

CWP 106H

Baxter

Gambro
CWP100
DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

The Gambro CWP 100 system is an integral part of the Hygienic Chain concept and allows for achieving and maintaining of the fluid purity³

- Automated heat disinfection of the distribution loop
- Fully enclosed central water plant
- Water saving functionality
- Chemical disinfection of the RO membranes
- Quiet, reliable operation

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CWP 106H

Baxter

Gambro
CWP100
DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

The Gambro CWP 100 system is an integral part of the Hygienic Chain concept and allows for achieving and maintaining of the fluid purity.³

- Automated heat disinfection of the distribution loop
- Fully enclosed central water plant
- Water saving functionality
- Chemical disinfection of the RO membranes
- Quiet, reliable operation

Making possible personal

Brochure

Download

» WATER SYSTEM

BACK TO OVERVIEW

WRO 300H

CWP 64 66

CWP 106H

Baxter

Gambro
CWP60

DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

- User-friendly control panel
- Simple, reliable and compact design
- Integration into Gambro dialysis machines
- Low noise level
- The CWP 60 system is available in five capacities³

Making possible personal.

The CWP 60 system

TECHNICAL DATA

PRODUCTION WATER

Model	64	66
at 5°C	23 (9120)	30 (11920)
+10°C	25 (9500)	32 (11920)
+15°C	27 (10200)	35 (11920)
+20°C	29 (11760)	35 (11920)

QUALITY

Dependent on inlet water quality. If pre-treated water is used and the system is properly maintained, the following rejection rates will be obtained:	
Total dissolved solids	> 95%
Bacteria (CFU) & endotoxins (EU)	> 99%
Water conversion factor	max. 75%

WATER SUPPLY

Model	64	66
Min. input, l/min (hour)	50 (2000)	60 (2400)
Min. input pressure, kPa	250	300
Max. input pressure	500 kPa (at model)	
Temperature	+5°C to +25°C	

QUALITY

Pre-treatment is recommended. Additional pre-treatment is normally necessary. Membrane and machine life expectancy depend on inlet water quality. It is recommended not to operate the CWP outside the following limits:	
Hardness	< 1.5 (6.0 °dH)
Iron	< 0.1 mg/l
Manganese	< 0.1 mg/l
Chloride	< 100 mg/l
Silica	< 20 mg/l
Total dissolved solids	< 1500 mg/l
Chlorine	< 0.1 mg/l
Fouling index (SIF 1)	< 5

CLEAN REQUIREMENTS

Model	64	66
Min. input, l/min (hour)	50 (2000)	60 (2400)

Maximum flow occurs during flushing.

CONNECTIONS

Inlet	PVC union, 8 to pipe OD 25 mm
Drain	PVC union, 8 to pipe OD 32 mm
Product water outlet & return	PVC union, 8 to pipe OD 25 mm

MEMBRANES

Membrane material	Modified polyamide, thin film composite
Membrane configuration	Spiral wound
pH tolerance	2-11

POWER SUPPLY

Main voltage	230/400 V, 50 Hz, three phase
Power rating	WRO 64, 66: 2.3 kW
Fuse	WRO 64, 66: 10 A (slow blow)

AMBIENT TEMPERATURE

Shipping and storage	+1°C to +41°C
Operation	+5°C to +31°C

DIMENSIONS

Model	64	66
Depth (mm)	410	410
Width (mm)	1030	1030
Height (mm)	2025	2025

80 mm extra space is required on both sides of the unit.

WEIGHT

Model	64	66
Unit (kg)	225	240
Unit + packing (kg)	360	390
Operation (kg)	240	300

1. Joseph M. et al. Blood Purif 2002;24:11-18. The Importance of Water Quality and Hemodialysis Fluid Composition.
2. ISO 15912:2015
3. User Manual CWP 60, 2012

CE 0066 This product is CE-marked in accordance with the requirements in EC Council Directive 93/42/EEC of 14 June 1993 concerning medical devices. For further information and operating instructions, please refer to applicable operator's manual.

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Baxter and Gambro are trademarks of Baxter International Inc. or its subsidiaries.

100000270000

» WATER SYSTEM

WRO 300H

Baxter

WRO300H

WATER PURITY FOR SINGLE PATIENTS

Automated procedures increase patient and caregiver involvement
 Dialysis machine can start and stop the Gambro WRO 300 H unit
 Automatic shut-off at completion of disinfection process
 The WRO 300 H unit can be programmed to automatically start heat disinfection at specified intervals


- Hot water disinfection
- Low noise level

Making possible personal.

Brochure

Download

CWP 64, 66

Baxter

Gambio

CWP60

DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}


- User-friendly control panel
- Simple, reliable and compact design
- Integration into Gambio dialysis machines
- Low noise level
- The CWP 60 system is available in five capacities³

Making possible personal.

Brochure

Download

CWP 106H

Baxter

Gambio

CWP100

DESIGNED FOR WATER PURITY

The quality of the water used to prepare dialysis fluid is extremely important. Even water considered as acceptable according to international standards may have chronic as well as acute effects on dialysis patients.^{1,2}

The Gambio CWP 100 system is an integral part of the Hygienic Chain concept and allows for achieving and maintaining of the fluid purity.³

- Automated heat disinfection of the distribution loop
- Fully enclosed central water plant
- Water saving functionality
- Chemical disinfection of the RO membranes
- Quiet, reliable operation

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REFERENCES [BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2/X? *Kidney International* 2015; 88: 480-485.
2. Terzani L, et al. CRRT guidance on dialysis strategies. *Neephrol Dial Transplant* 2007; 22 Suppl_2: i5-i25.
3. Terzani L, et al. Online HDF definition, dose quantification and safety revisited. *Neephrol Dial Transplant* 2012; 27: S42-S45.
4. Wiley M et al. Examining the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3): 2213-2219; quiz 2013(975).
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Curr Opin Nephrol Dial* 2017; 19: 142-152.
6. Rocco C. The role of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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* Do not use **TheraNova** dialyzers in HDF or HF mode

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REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2V? *Kidney International* 2015; 88: 480-485.
2. Teramachi J, et al. CRRT guidance on dialysis strategies. *Neophyl Dial Transplant* 2007; 22 Suppl 2: 15-125.
3. Teramachi J, et al. Online HDF definition, dose quantification and safety revisited. *Neophyl Dial Transplant* 2012; 27: 542-555.
4. Wiley M et al. Examining the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3):2213-2219.
5. Hutchman CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Contrib Nephrol* 2017; 193:142-152.
6. Rocco C. The rise of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2V? *Kidney International* 2020; 98: 480-485.
2. Teramachi J, et al. CRRT guidance on dialysis strategies. *Neonatal Dial Transplant* 2007; 22 Suppl 2: 45-125.
3. Teramachi J, et al. Online HDF definition, dose quantification and safety revisited. *Neonatal Dial Transplant* 2012; 27: 542-555.
4. Wilby M et al. Efficacy of the Clinical Efficacy of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3): 2213-2219; quiz 2016; 10: 1365-71.
5. Hutchman CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Crit Care Med* 2017; 45: 142-152.
6. Rocco C. The rise of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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[Read more](#)

* Do not use TheraNova dialyzers in HDF or HF mode

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* Do not use TheraNova dialyzers in HDF or HF mode

REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2/Y? *Kidney International* 2015; 88: 480-485.
2. Terrence J, et al. Online Kt/V: a dialysis adequacy measure. *Hepathol Dial Transplant* 2007; 22 Suppl 2: 45-52.
3. Terrence J, et al. Online Kt/V: definition, dose quantification and safety revisited. *Hepathol Dial Transplant* 2012; 22: 542-555.
4. Wible M et al. Examining the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3):2221N/CJN.10130515.
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Contrib Nephrol* 2017; 193:142-152.
6. Rocco C. The rise of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

MONITOR

- AK 98 system

TREATMENT MONITORING

- Diascan
- BPM

WATER SYSTEM

- CWP 100 with integrated heat disinfection

TECHNICAL SERVICES

- Hotline
- Technical training
- Service contract

* Do not use TheraNova dialyzers in HDF or HF mode

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OVERALL CLEARANCE HD vs. HDx vs. HDF

HD with latest generation high-flux dialyzer


HDx with Theranova 400 dialyzer

HDF with latest generation high-flux dialyzer for HDF

* p<0.001 vs Theranova 400

** p<0.01 vs Theranova 400

Q_b = 400 ml/min - Treatment Time = 4.4 h - Vconv = 24L (Mean) - n = 20

Adapted from 14) Kirsch et al. *Nephrol Dial Transplant* 2017

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REFERENCES

BACK TO OVERVIEW

1. Venkatesh R, et al. Once upon a time in dialysis: the last days of RRT. *Kidney International* 2015; 88: 480-485.
2. Teramachi J, et al. Online CRRT guidance on dialysis strategies. *Neperol Dial Transplant* 2007; 22 Suppl 2: 45-025.
3. Teramachi J, et al. Online HDQ definition, dose quantification and safety revisited. *Neperol Dial Transplant* 2012; 27: 542-055.
4. Wille M et al. Examining the Clinical Relevance of Predicting Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3): 10.2215/CJN.10131615.
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Curr Opin Nephrol Dialysis* 2017; 19: 142-152.
6. Rocco C. The Use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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REFERENCES [BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of RFX? *Kidney International* 2015; 88: 480-485.
2. Terresato J, et al. CRRT guidance on dialysis strategies. *Intensive Care Medicine* 2017; 22 Suppl 2: 45-52.
3. Terresato J, et al. Online HDF definition, dose quantification and safety revisited. *Intensive Care Medicine* 2012; 22: 542-555.
4. Wiley M et al. Extending the Clinical Evidence of Proving Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3): 221N/CN.10113615.
5. Hutchison CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Crit Care Med* 2017; 45: 142-152.
6. Rocco C. The role of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2/Y? *Kidney International* 2015; 88: 480-485.
2. Terzani J, et al. CRRT guidance on dialysis strategies. *Neophrol Dial Transplant* 2007; 22 Suppl_2: 45-125.
3. Terzani J, et al. Online HDF definition, dose quantification and safety revisited. *Neophrol Dial Transplant* 2012; 27: 542-555.
4. Wille M et al. Efficacy of the Clinical Efficacy of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3):2215-2221; quiz 2015(3).
5. Hutchman CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Crit Care Med* 2017; 45:142-152.
6. Rocco C. The rise of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2/KT. *Kidney International* 2015; 88: 480-485.
2. Terrence J, et al. CRRT guidance on dialysis strategies. *Intensive Care Medicine* 2017; 22 Suppl_2: 45-52.
3. Terrence J, et al. Online HDF definition, dose quantification and safety revisited. *Intensive Care Medicine* 2012; 22: 542-555.
4. Wiley M et al. Examining the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol* 2016; Mar 9 (3):2213/CJN.10130515.
5. Hutchison CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Contrib Nephrol* 2017; 193:142-152.
6. Rocco C. The Use of Expanded Hemodialysis. *Blood Purif* 2017; 44: L-VI.

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[Read more](#)

* Do not use Theraflow dialyzers in HDF or HF mode.

Baxter

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Provide Expanded HD. Retain HD simplicity

* Do not use Theraflow dialyzers in HDF or HF mode.

REFERENCES

[BACK TO OVERVIEW](#)

1. Venkater S, et al. Once upon a time in dialysis: the last days of R2V? *Kidney International* 2015; 88: 480-485.
2. Terzani L, et al. CRP, guanine on dialysis strategies. *Neperol Dial Transplant* 2007; 22 Suppl 2: 45-025.
3. Terzani L, et al. Online HDF definition, dose quantification and safety revisited. *Neperol Dial Transplant* 2012; 22: 542-055.
4. Wille M et al. Extending the Clinical Evidence of Providing Increased Removal of Large Middle Molecules. *Dis J Am Soc Nephrol* 2016; Mar 9 Issue 10:2215/CLIN.10130515.
5. Hutchison CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Curr Opin Nephrol* 2017; 19:142-152.
6. Rocco C. The role of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-16.

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REFERENCES

[BACK TO OVERVIEW](#)

7. Bunch A et al. *Long Term Effects of Expanded Hemodialysis (HDx) on Clinical and Laboratory Parameters in a Large Cohort of Dialysis Patients*. ASN 2018 Kidney Week: FR-PO766.
8. Chapdelaine I, et al. *Optimization of the convection volume in online post-dilution haemodiafiltration: practical and technical issues*. Clinical Kidney Journal 2015; 8: 191-198.
9. Hutchison CA, et al. *Quantitative assessment of serum and urinary polyclonal free light chains in patients with chronic kidney disease*. Clin J Am Soc Nephrol 2008; 3:1684-90.
10. Duranton F, et al. *Normal and pathologic concentrations of uremic toxins*. J Am Soc Nephrol 2012; 23(7):1258-70.
11. Okyay GU, et al. *Novel inflammatory marker in dialysis patients: YKL-40*. Ther Apher Dial 2013; 17: 193-201.
12. Cobo G, et al. *Chronic inflammation in end-stage renal disease and dialysis*. Nephrol Dial Transplant 2018 33: iii35–iii40


[Back to Portfolio Vablet](#)

REFERENCES

[BACK TO OVERVIEW](#)

13. Sakurai K. *Biomarkers for evaluation of clinical outcomes of hemodiafiltration*. Blood Purif 2013; 35 Suppl_1: 64-8.
14. Kirsch AH, et al. *Performance of hemodialysis with novel medium cut-off dialyzers*. Nephrol Dial Transplant 2017; 32: 165-172.
15. Ronco C, et al. *Evolution of synthetic membranes for blood purification: the case of the Polyflux family*. Nephrol Dial Transplant (2003) 18 Suppl_7: vii10–vii20
16. Schepers E, et al. *Assessment of the association between increasing membrane pore size and endotoxin permeability using a novel experimental dialysis simulation set-up*. BMC Nephrology 2018; 19:1.
17. Baxter Extract from Operator Manual AK98_2xx_HCEN12745 10/2016 pp A24 & A91-94.
18. Data on file "*Therapy switching cost*" from *Baxter HDx Cost calculator*" EUMP/MG209/16-0002a.


[Back to Portfolio Vablet](#)

REFERENCES

[BACK TO OVERVIEW](#)

19. Data on file "Therapy cost per treatment" from Baxter HDx Cost calculator EUMP/MG209/16-0002a

31. Yilmaz MI, et al. *Low-grade inflammation in chronic kidney disease patients before the start of renal replacement therapy: sources and consequences* Clin Nephrol 2007; 68: 1-9.

32. Stenvinkel P. *Can Treating Persistent Inflammation Limit Protein Energy Wasting?* Seminars in Dialysis 2013; 26(1): 16–39.

33. Akchurin OM, Kaskel F. *Update on Inflammation in Chronic Kidney Disease.* Blood Purif 2015; 39: 84–92.

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Baxter, AK 98, BiCart, Diascan, CleanCart, Novaline and Theranova are trademarks of Baxter International Inc. or its subsidiaries.

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* Do not use TheraSense dialysis in HCF or HF mode

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* Do not use TheraSense dialysis in HCF or HF mode

Baxter

REFERENCES

[BACK TO OVERVIEW](#)

1. Kishimoto R, et al. Once upon a time in dialysis: the last days of BIC? *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESPD guideline on dialysis strategies. *Neperol Dial Transplant* 2007; 22 Suppl 2: i6-i21.
3. Tattersall J, et al. Drive HD? dialysis, dose quantification and safety revisited. *Neperol Dial Transplant* 2012; 22: 542-550.
4. Wiley M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule. *Clin J Am Soc Nephrol* 2018; Mar 8 (vol 12): 215-221; quiz 216-217.
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HEM). *Clinical Nephrol* 2017; 193:142-150.
6. Bello C. The Use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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* Do not use Theranova dialysis in HCF or HF mode

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* Do not use Theranova dialysis in HCF or HF mode

REFERENCES

[BACK TO OVERVIEW](#)

1. Karkhanavala R, et al. Once upon a time in dialysis: the last days of RDT? *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESRD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl 2: i6-i21.
3. Tattersall J, et al. Online HDF: definition, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 542-550.
4. Wible M et al. Expanding the Clinical Reference of Providing Increased Removal of Large Molecule Wastes. *Clin J Am Soc Nephrol* 2018; Mar 8 (vol 12): 215-221; quiz 216-217.
5. Hutchinson CA, et al. The Advantage for Expanded Hemodialysis Therapy (HEM). *Clinical Nephrology* 2017; 193:142-150.
6. Benic C. The New of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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
Adapted from 18) Caplin et al. *Nephrol Dial Transplant*, 2011

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[Read more](#) [Read more](#) [Read more](#)

* Do not use Theranova dialyzers in HDF or HF mode

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CONSUMABLES

SERVICES & SYSTEM

* Do not use Theranova dialyzers in HDF or HF mode

REFERENCES

BACK TO OVERVIEW

1. Kishore R, et al. Once upon a time in dialysis: the last days of BIC? *Kidney International* 2015; 88: 460-465.
2. Tavesani J, et al. EPO guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl 2: i6-i21.
3. Tavesani J, et al. Online HDF: definition, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 942-950.
4. Wille M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule Wastes. *Clin Am Soc Nephrol* 2018; Mar 5 (vol 10): 2215-2216. DOI: 10.1093/cn.109.75
5. Hutchinson CA, et al. The Advantage for Expanded Hemodialysis Therapy (HEM). *Clinical Nephrology* 2017; 193:142-150.
6. Senz C. The use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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
* Do not use Theraflow dialyzers in HDF or HF mode

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MONITOR: • Artis Physio Plus

TREATMENT MONITORING: • Diascan • BPM • HemoControl

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* Do not use Theraflow dialyzers in HDF or HF mode

Baxter

REFERENCES

BACK TO OVERVIEW

1. Santoro R, et al. Once upon a time in dialysis: the last days of RDT. *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESPD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl 2: i6-i21.
3. Tattersall J, et al. Online HDF: definition, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 542-550.
4. Wille M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule. *Clin J Am Soc Nephrol* 2018; Mar 8 (vol 10): 2215-2216, DOI:10.1093/cjn.
5. Hutchinson CA, et al. The Advantage for Expanded Hemodialysis Therapy (HEM). *Clinical Nephrology* 2017; 193:142-150.
6. Bello C. The use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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CONSUMABLES MONITOR & FEEDBACK SERVICES & SYSTEM

HDx-HC

REFERENCES

BACK TO OVERVIEW

1. Karkhanavala R, et al. Once upon a time in dialysis: the last days of HD? *Kidney International* 2015; 88: 460-465.
2. Tavesani J, et al. ESRD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl_2: i6-i21.
3. Tavesani J, et al. Online HD? adhesion, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 942-950.
4. Wille M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule. *Clin J Am Soc Nephrol* 2018; Mar 5 (vol 10): 2215-2216; quiz 2216-2217.
5. Hutchinson CA, et al. The Advantage for Expanded Hemodialysis Therapy (HEM). *Current Nephrol* 2017; 19:1342-1350.
6. Somo G. The use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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* Do not use Theranova/dialysate in HCF or HF mode

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* Do not use Theranova/dialysate in HCF or HF mode

REFERENCES

[BACK TO OVERVIEW](#)

1. Kishimoto R, et al. Once upon a time in dialysis: the last days of BIC®. *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESPD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl 2: ii-421.
3. Tattersall J, et al. Drive HD® adhesion, dose quantification and safety verified. *Kidney Dial Transplant* 2012; 27: 942-950.
4. Wible M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule Wastes. *Clin J Am Soc Nephrol* 2018; Mar 8 (vol 12): 215-CAK.12018.75
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HEM). *Clinical Nephrol* 2017; 193:142-150.
6. Somo G. The use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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
REFERENCES

[BACK TO OVERVIEW](#)

1. Kurokawa R, et al. Once upon a time in dialysis: the last days of RDT? *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESPD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl_2: i6-i21.
3. Tattersall J, et al. Online HDF adhesion, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 542-550.
4. Willeke M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule. *Clin J Am Soc Nephrol* 2018; Mar 8 (vol 10): 2215; DOI:10.1093/cjn.
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HEP). *Clinical Nephrology* 2017; 193:142-150.
6. Somo G. The use of Expanded Hemodialysis. *Blood Purif* 2017; 44:1-VII.

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* Do not use TheraNova-dialyzer in HCF or HF mode

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SERVICES & SYSTEM

- HDx-HC
- HDx-HC
- HDx-HC

* Do not use TheraNova-dialyzer in HCF or HF mode

REFERENCES

[BACK TO OVERVIEW](#)

1. Kishimoto R, et al. Dose upon a time in dialysis: the last days of Kt/V? *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESPD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl 2: i6-i21.
3. Tattersall J, et al. Online HDF: definition, dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 142-150.
4. Willeke M et al. Expanding the Clinical Relevance of Proving Increased Removal of Large Molecule. *Clin Am Soc Nephrol* 2018; Mar 5 (vol 10): 2215-2216, CAK.L10.100675
5. Hutchinson CA, et al. The Advantage for Expanded Hemodialysis Therapy (HDx). *Clinical Nephrology* 2017; 193:142-150.
6. Ferris C. The Use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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[Read more](#)

* Do not use Theranova® dialysis in HCF or HF mode

Baxter

THE HDx-HC TREATMENT OPTION

Provide expanded HD. Reduce intradialytic hypotension

* Do not use Theranova® dialysis in HCF or HF mode

REFERENCES

[BACK TO OVERVIEW](#)

1. Karkhanavich R, et al. Once upon a time in dialysis: the last days of BIC®. *Kidney International* 2015; 88: 460-465.
2. Tattersall J, et al. ESRD guideline on dialysis strategies. *Kidney Dial Transplant* 2007; 22 Suppl_2: i6-i21.
3. Tattersall J, et al. Drive HD® dialysis: dose quantification and safety revisited. *Kidney Dial Transplant* 2012; 27: 942-950.
4. Wiley M et al. Expanding the Clinical Relevance of Providing Increased Removal of Large Molecule Toxins. *Clin J Am Soc Nephrol* 2018; Mar 5 (suppl 3): S215-GRA.1311.09675.
5. Hutchinson CA, et al. The Rationale for Expanded Hemodialysis Therapy (HDx). *Current Nephrol* 2017; 19(1):142-150.
6. Somo G. The Use of Expanded Hemodialysis. *Blood Purif* 2017; 44: 1-10.

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REFERENCES

[BACK TO OVERVIEW](#)

1. Vanholder R, et al. *Once upon a time in dialysis: the last days of Kt/V?*. *Kidney International* 2015; 88: 460–465.
2. Tattersall J, et al. *EBPG guideline on dialysis strategies*. *Nephrol Dial Transplant* 2007; 22 Suppl_2: ii5–ii21.
3. Tattersall J, et al. *Online HDF definition, dose quantification and safety revisited*. *Nephrol Dial Transplant* 2012; 22: 542–550.
4. Wolley M et al. *Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules*. *Clin J Am Soc Nephrol*. 2018; Mar 5 (doi:10.2215/CJN.10110917).
5. Hutchison CA, et al. *The Rationale for Expanded Hemodialysis Therapy (HDx)*. *Contrib Nephrol* 2017; 191:142-152.
6. Ronco C. *The rise of Expanded Hemodialysis*. *Blood Purif* 2017; 44: I–VIII.

[Back to Portfolio Vablet](#)

REFERENCES

[BACK TO OVERVIEW](#)

7. Bunch A et al. *Long Term Effects of Expanded Hemodialysis (HDx) on Clinical and Laboratory Parameters in a Large Cohort of Dialysis Patients*. ASN 2018 Kidney Week: FR-P0766.
8. Chapdelaine I, et al. *Optimization of the convection volume in online post-dilution haemodiafiltration: practical and technical issues*. Clinical Kidney Journal 2015; 8: 191-198.
9. Hutchison CA, et al. *Quantitative assessment of serum and urinary polyclonal free light chains in patients with chronic kidney disease*. Clin J Am Soc Nephrol 2008; 3:1684-90.
10. Duranton F, et al. *Normal and pathologic concentrations of uremic toxins*. J Am Soc Nephrol 2012; 23(7):1258-70.
11. Okyay GU, et al. *Novel inflammatory marker in dialysis patients: YKL-40*. Ther Apher Dial 2013; 17: 193-201.
12. Cobo G, et al. *Chronic inflammation in end-stage renal disease and dialysis*. Nephrol Dial Transplant 2018 33: iii35–iii40

[Back to Portfolio Vablet](#)

REFERENCES

[BACK TO OVERVIEW](#)

13. Sakurai K. *Biomarkers for evaluation of clinical outcomes of hemodiafiltration*. Blood Purif 2013; 35 Suppl_1: 64-8.
14. Kirsch AH, et al. *Performance of hemodialysis with novel medium cut-off dialyzers*. Nephrol Dial Transplant 2017; 32: 165-172.
15. Ronco C, et al. *Evolution of synthetic membranes for blood purification: the case of the Polyflux family*. Nephrol Dial Transplant (2003) 18 Suppl_7: vii10–vii20
16. Schepers E, et al. *Assessment of the association between increasing membrane pore size and endotoxin permeability using a novel experimental dialysis simulation set-up*. BMC Nephrology 2018; 19:1.
17. Baxter Extract from Operator Manual Artis Physio Plus 2018 code 9033282100 pages 175-182.
18. Caplin B, et al. *Patients' perspective of haemodialysis-associated symptoms*. Nephrol Dial Transplant 2011; 26(8): 2656-2663.

Back to Portfolio Vablet

REFERENCES

BACK TO OVERVIEW

19. Stefánsson B, et al. *Intradialytic Hypotension and Risk of Cardiovascular Disease*. Clin J Am Soc Nephrol 2014; 9: 2124–2132.
20. Winkler RE, et al. *Blood Volume Regulation*. In: Technical Problems in Patients on Hemodialysis. Rijeka. Croatia. 2011: 235-250.
21. *National Kidney Foundation Kidney Disease Outcomes Quality Initiative. Clinical practice guidelines for cardiovascular disease in dialysis patients*. Available at: http://www2.kidney.org/professionals/KDOQI/guidelines_cvd/intradialytic.htm. Accessed: April 2017.
22. Selby N, et al. *Acute cardiac effect of dialysis*. Seminars in Dialysis 2007; 20(3): 220–228.
23. Gayraud N, et al. *Consequences of increasing convection onto patient care and protein removal in hemodialysis*. PLoS ONE 2017; 12(2): e0171179.
24. Kooman J, et al. *EBPG guideline on haemodynamic instability*. Nephrol Dial Transplant 2007; 22 Suppl_2: ii22–ii44.

[Back to Portfolio Vablet](#)

REFERENCES

[BACK TO OVERVIEW](#)

- 25. Nesrallah GE, et al. *Biofeedback dialysis for hypotension and hypervolemia: a systematic review and meta-analysis.* Nephrol Dial Transplant. 2013; 28:182–191.

- 26. Basile C, et al. *Efficacy and safety of haemodialysis treatment with the Hemocontrol biofeedback system: a prospective medium-term study.* Nephrol Dial Transplant. 2001; 16(2): 328-34.

- 27. Doria M, et al.. *The dialysis staff workload and the blood volume tracking system during the hemodialysis sessions of hypotension-prone patients* Int J Artif Organs. 2014; 37(4): 292–298.

- 28. Gil HW, et al. *Efficacy of hemocontrol biofeedback system in intradialytic hypotension-prone hemodialysis patients.* J Korean Med Sci. 2014; 29: 805–810.

- 29. Deziel C, et al. *Impact of hemocontrol on hypertension, nursing interventions, and quality of life: A randomised, controlled trial.* Clin J Am Soc Nephrol. 207; 2: 661–668.

- 30. Selby N, et al. *Occurrence of Regional Left Ventricular Dysfunction in Patients Undergoing Standard and Biofeedback Dialysis.* Am J Kidney Dis 2006; 47: 830-841.

[Back to Portfolio Vablet](#)

[BACK TO OVERVIEW](#)

REFERENCES

- 31. Yilmaz MI, et al. *Low-grade inflammation in chronic kidney disease patients before the start of renal replacement therapy: sources and consequences.* Clin Nephrol 2007; 68: 1-9.

- 32. Stenvinkel P. *Can Treating Persistent Inflammation Limit Protein Energy Wasting?* Seminars in Dialysis 2013; 26(1): 16-39.

- 33. Akchurin OM, Kaskel F. *Update on Inflammation in Chronic Kidney Disease.* Blood Purif 2015; 39: 84-92.

For the safe and proper use of the devices mentioned herein, refer to the Instructions For Use or Operator's Manual.

Baxter, Artis Physio, Artiset, Diascan, HemoControl, CleanCart, and Theranova are trademarks of Baxter International Inc. or its subsidiaries.