Baxter

HDx

Enabled by
MCO THERANOVA

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Health resource utilization	THERANOVA [n = 86]	high-flux HD [n = 85] ^a
Hospitalization events	18	31
Total hospital days	74	139
Total patient-years	32.4	30.5
Hospitalization rate per PY [SE]	0.56 [0.13]	1.02 [0.12]
Hospital length of stay [mean days [SE]]	4.11 [0.57]	4.63 [0.58]

^a One high-flux HD randomized participant did not complete baseline.

THERANOVA dialyzers are indicated for treatment of chronic and acute renal failure by Hemodialysis

Do not use in Hemodiafiltration or Hemofiltration mode or isolated Ultrafiltration

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
Uremic Toxin Classes by molecular weight [Daltons]^{32,50}

Urea	[60 Da]	●	Small Molecules [<0.5 kDa]
Phosphate	[96 Da]	●	
PTH	[9.5 kDa]	●	
Beta2 microglobulin	[12 kDa]	●	
Myoglobin	[17 kDa]	●	
Kappa free-light-chains	[23 kDa]	●	
Complement factor D	[24 kDa]	●	
Interleukin-6	[25 kDa]	●	
TNF-alpha	[26 kDa]	●	
FGF-23	[32 kDa]	●	
Alpha 1 microglobulin	[33 kDa]	●	
YKL-40	[40 kDa]	●	
Lambda free-light-chains	[45 kDa]	●	
Albumin	[67 kDa]	●	

Evolution of dialysis therapies

Adapted after Rosner M, et al. Classification of Uremic Toxins and Their Role in Kidney Failure.

Clin J Am Soc Nephrol. 2021;16(12):1918-1928³²

EUTOX Uremic Solutes Database. June 2022. Uremic-toxins.org⁵⁰

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
Large-middle molecule linkage with clinical symptoms and outcomes

Large-middle molecules		Relevant clinical effects
TNF-alpha	[26 kDa]	<ul style="list-style-type: none"> Sepsis³³ Chronic Inflammation³³ Cardiovascular Disease³⁴ Protein-energy wasting in CKD³⁴
FGF-23 ⁵⁰	[32 kDa]	<ul style="list-style-type: none"> Secondary Immunodeficiency Cardiovascular Disease³⁴
Alpha 1 microglobulin	[33 kDa]	<ul style="list-style-type: none"> Restless Legs Syndrome [RLS]^{35,36}
YKL-40	[40 kDa]	<ul style="list-style-type: none"> Inflammation³⁷
Lambda free-light-chains	[45 kDa]	<ul style="list-style-type: none"> Chronic Inflammation Secondary Immunodeficiency³⁴

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Number of pores

Adapted after Wolley: Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules.⁴⁹

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
Adapted after Boschetti-de-Fierro: MCO Membranes: Enhanced Selectivity in High-Flux Class.³⁸

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
THERANOVA dialyzers are indicated for treatment of chronic and acute renal failure by Hemodialysis.
Do not use in Hemodiafiltration or Hemofiltration mode or isolated Ultrafiltration
Medical devices of class IIb – Notified body: BSI, NL [CE 2797] – Legal manufacturer: Gambro Dialysatoren GmbH – Hechingen, Germany.¹
For single use only. For safe and proper use of these devices refer to the Instructions for Use.

REFERENCES

1. Zweigart C, Boschetti-de-Fierro A, Hulko H, et al. Medium cut-off membranes – closer to the natural kidney removal function. *Int J Artif Organs.* 2017;40(7):328-334.
2. Ronco, C. The Rise of Expanded Hemodialysis. *Blood Purif.* 2017;44:I-VIII.
3. Molano-Trivino A, Sanabria M, Vesga J, Buitrago G, Sanchez R, Rivera A. Effectiveness of medium cut-off vs high flux dialyzers: a propensity score matching cohort study. *In Nephrol Dial Transport.* 2021;36:486-U948.
4. Penny JD, Jarosz P, Salerno FR, Lemoine S, McIntrye CW. Impact of Expanded Hemodialysis Using Medium Cut-off Dialyzer on Quality of Life: Application of Dynamic Patient-Reported Outcome Measurement Tool. *Kidney Medicine.* 2021;3(6):992-1002.
5. Hutchison CA, Wolley M. The Rationale for Expanded Hemodialysis Therapy [HDx]. *Contrib Nephrol.* 2017;191:142-52.
6. Kirsch AH, Lyko M, Nilsson LG, et al. Performance of hemodialysis with novel medium cut-off dialyzers. *Nephrol Dial Transpl.* 2017;32(1):165-72.
7. Zhang JC, El-Majzoub S, Li M, et al. Could symptom burden predict subsequent healthcare use in patients with end stage kidney disease on hemodialysis care? A prospective, preliminary study. *Ren Fail.* 2020;42(1):294-301.
8. Lim JH, Park Y, Yook JM, et al. Randomized controlled trial of medium cut-off versus high-flux dialyzers on quality of life outcomes in maintenance hemodialysis patients. *Sci Rep.* 2020;10(1):1-11.
9. Alarcon J.C, Bunch A, Ardila F, et al. Impact of Medium Cut-Off Dialyzers on Patient-Reported Outcomes: COREXH Registry. *Blood Purif.* 2021;50:110-118.
10. Anwar N, Naz N, Reynolds A, Chamber. HDx: is it a better way to dialyse? *Nephrol Dial Transplant.* 2020;35(3):1109.
11. Bolton S, Gair R, Nilsson LG, Matthews M, Stewart L, McCullagh N. Clinical Assessment of Dialysis Recovery Time and Symptom Burden: Impact of Switching Hemodialysis Therapy Mode. *Patient Relat Outcome Meas.* 2021;12:315-321.
12. Rayner HC, Zepel L, Fuller DS, et al. Recovery time, quality of life, and mortality in hemodialysis patients: the Dialysis Outcomes and Practice Patterns Study [DOPPS]. *Am J Kidney Dis.* 2014;64:86-94.
13. Cozzolino M, Magagnoli L, Ciceri P, Conte F, Galassi A. Effects of a medium cut-off [THERANOVA®] dialyser on haemodialysis patients: a prospective, cross-over study. *Clin. Kidney J.* 2021;14(1):382-389.
14. Sanabria RM, Hutchison CA, Vesga, JL, Ariza JG, Sanchez R, Suarez AM. Expanded Hemodialysis and Its Effects on Hospitalizations and Medication Usage: A Cohort Study. *Nephron.* 2021;145(2):179-187.
15. Lim JH, Jeon Y, Yook JM, et al. Medium cut-off dialyzer improves erythropoiesis stimulating agent resistance in a hepcidin-independent manner in maintenance hemodialysis patients: results from a randomized controlled trial. *Sci Rep.* 2020;10(1):1-10.
16. Ariza JG, Walton SM, Suarez AM, Sanabria M, Vesga JL. An initial evaluation of expanded hemodialysis on hospitalizations, drug utilization, costs, and patient utility in Colombia. *Ther Apher Dial.* 2021;25(5):621-627.
17. Hadad-Arrascue F, Nilsson LG, Rivera AS, Bernardo AA, Cabezaudo Romero JB. Expanded hemodialysis as effective alternative to on-line hemodiafiltration: A randomized mid-term clinical trial. *Ther Apher Dial.* 2022;26(1):37-44.
18. Pisoni RL, Wikstrom B, Elder SJ, et al. Pruritus in haemodialysis patients: International results from the Dialysis Outcomes and Practice Patterns Study [DOPPS]. *Nephrol Dial Transplant.* 2006;21(12):3495-3505.
19. Su H, Lei C.T, Zhang C. Interleukin-6 Signaling Pathway and Its Role in Kidney Disease: An Update. *Front Immunol.* 2017;8:405.
20. Kimmel, M, Alischer OM, Dunst R, Braun, et al. The role of micro-inflammation in the pathogenesis of uraemic pruritus in haemodialysis patients. *Nephrology Dialysis Transplantation.* 2006;21(3):749-755.
21. Lin XW, Zhang JF, Qiu MY, et al. Restless legs syndrome in end stage renal disease patients undergoing hemodialysis. *BMC Neurol.* 2019;19:47.
22. Freire de Menezes A, et al. *Int J Nephrol.* 2018;25:2018:1414568.
23. Giannaki CD, et al. Epidemiology, impact and treatment options of restless legs syndrome in end-stage renal disease patients; an evidence based review. *Kidney Int.* 2014;85(6):1275-1282.
24. La Manna G, Pizza F, Persici E, et al. Restless legs syndrome enhances cardiovascular risk and mortality in patients with end-stage kidney disease undergoing long-term haemodialysis treatment. *Nephrol Dial Transplant.* 2011;26:1976-1983.
25. Olsson M.G, Allhorn M, Bulow L, et al. Pathological conditions involving extracellular hemoglobin molecular mechanisms, clinical significance, and novel therapeutic opportunities for alpha-1-microglobulin. *Antioxid Redox Signal.* 2012;17:813-846.
26. Bossola M, Tazza L. Postdialysis fatigue: a frequent and debilitating symptom. *Semin Dial.* 2016;29:222-227.
27. Bossola M, Vulpio C, Tazza L. Fatigue in chronic dialysis patients. *Semin Dial.* 2011;24:550-555.
28. Lindsay RM, Heidenheim PA, Nesrallah G, Garg AX, Suri R. Daily hemodialysis study group London health sciences centre: Minutes to recovery after a hemodialysis session: a simple health-related quality of life question that is reliable, valid, and sensitive to change. *Clin J Am Soc Nephrol.* 2006;1:952-959.
29. Bossola M, Di Stasio E, Monteburini T, et al. Recovery Time after Hemodialysis Is Inversely Associated with the Ultrafiltration Rate. *Blood Purif.* 2019;47(1-3):45-51.
30. Bossola M, Di Stasio E, Giungi S, Rosa F, Tazza L. Fatigue is associated with serum interleukin-6 levels and symptoms of depression in patients on chronic hemodialysis. *J Pain Symptom Manag.* 2015;49:578-585.
31. Blackowicz MJ, Falzon L, Beck W, Tran H, Weiner DE. Economic evaluation of expanded hemodialysis with the THERANOVA 400 dialyzer: A post hoc evaluation of a randomized clinical trial in the United States. *Hemodialysis International.* 2022. <https://doi.org/10.1111/hdi.13015>.
32. Rosner M, Reis T, Husain-Syed, et al. Classification of Uremic Toxins and Their Role in Kidney Failure. *Clin J Am Soc Nephrol.* 2021;16(12):1918-1928.
33. Ronco C, Marchionna N, Brendolan A, Neri M, Lorenzin A, Martinez Rueda AJ. Expanded haemodialysis: from operational mechanism to clinical results. *Nephrology Dialysis Transplantation.* 2018;33(3):iii41-iii47.
34. Wolley M, Jardine M, Hutchison CA. Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol.* 2018;13(5):805-814.
35. Ronco, C. Expanded Hemodialysis: Innovative Clinical Approach in Dialysis. Vol.191. Karger Medical and Scientific Publishers. 2017.
36. Sakurai K. Biomarkers for Evaluation of Clinical Outcomes of Hemodiafiltration. *Blood Purif.* 2013;35(1):64-68.
37. Lorenz G, Schmalenberg M, Kemmner S, et al. Mortality prediction in stable hemodialysis patients is refined by YKL-40, a 40-kDa glycoprotein associated with inflammation. *Kidney Int.* 2018;93(1):221-230.
38. Boschetti-de-Fierro A, Voigt M, Storr M, et al. MCO Membranes: Enhanced Selectivity in High-Flux Class. *Sci Rep.* 2015;5(1):18448.
39. Ronco C, Clark WR. Haemodialysis membranes. *Nat Rev Nephrol.* 2018;14(6):394-410.
40. Device Classification Under Section 513(f)(2) [De Novo]. Accessdata.fda.gov. <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/de novo.cfm?id=DEN190042>. Published 2020. Accessed April 26, 2022.
41. Mazairac AH, Blankestijn PJ, Grootenhuis MP, et al. The cost-utility of haemodiafiltration versus haemodialysis in the Convective Transport Study. *Nephrol Dial Transplant.* 2013;28(7):1865-1873.
42. Schepers E, Glorieux G, Eloot S, et al. Assessment of the Association Between Increasing Membrane Pore Size and Endotoxin Permeability Using a Novel Experimental Dialysis Simulation Set-Up. *BMC Nephrol.* 2018;19:1.
43. Boschetti-de-Fierro A, Beck W, Krause B, Hildebrandt H. Membrane Innovation in Dialysis. *Contrib Nephrol.* 2017;191:100-114.
44. Lorenzin A, Neri M, Lupi A, et al. Quantification of internal filtration in hollow fiber hemodialyzers with medium cut-off membrane. *Blood Purification.* 2018;46(3):196-204.
45. Zuo L, Wu J, Yao Q, et al. Classification of hemodialyzer clinical performance. *Zhongguancun Nephropathy Blood Purification Innovation Alliance.* 2021;1-4.
46. Clinical Data. Baxter Renal Care. <https://renalcare.baxter.com/data-insights/clinical-data>. Accessed April 26, 2022.
47. Kandi M, Brignardello-Petersen R, Couban R, Wu C, Nesrallah G. Clinical Outcomes with Medium Cut-Off Versus High-Flux Hemodialysis Membranes: A Systematic Review and Meta-Analysis. *Can J Kidney Health Dis.* 2022;9:1-16. DOI: 10.1177/20543581211067087
48. Ronco C, La Manna G. Expanded Hemodialysis: A New Therapy for a New Class of Membranes. *Contrib Nephrol.* 2017;190:124-133. doi:10.1159/000468959
49. Wolley M, Jardine M, Hutchison CA. Exploring the Clinical Relevance of Providing Increased Removal of Large Middle Molecules. *Clin J Am Soc Nephrol.* 2018;13(5):805-814. doi:10.2215/CJN.10110917
50. EUTOX Uremic Solutes Database. June 2022. [Uremic-toxins.org](http://uremic-toxins.org)
51. Wolley M, et al. *Clin J Am Soc Nephrol.* 2018;b;13(5):805-814.