

Baxter

Novaline

BLOODLINE TUBING SYSTEM

BLOODLINE TUBING SYSTEM

SIMPLE

EFFICIENT

NOVALINE MODELS FOR HD ON AK SYSTEMS

CODE	MODEL	SYSTEM	DN/SN	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE- DIALYZER CHAMBER	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955419	BL 05R	AK 95 S, AK 96, AK 98	DN	24	Ebeam	DEHP Free	8x12 mm	Yes	No	1	Yes	1	Dual	No
955312	BL 90 D4	AK 95 S, AK 96, AK 98	DN	24	Ebeam	DEHP Free	8x12 mm	No	No	1	Yes	1	Dual	No

NOVALINE MODELS FOR HD, HF AND HDF ON AK SYSTEMS

CODE	MODEL	SYSTEM	QUANTITY PER BOX	STERILIZATION	MATERIAL	PUMP SEGMENT	INFUSION LINE	BVS + COMBI CONNECTOR	TRANSDUCER PROTECTOR	PRE- DIALYZER CHAMBER	PRE- DILUTION	POST- DILUTION	VENOUS SERVICE LINE	DIALYZER CONNECTOR	REMOVABLE LABELING
955313	BL 90 D4T	AK 200 S, AK 200 Ultra S	24	Ebeam	DEHP Free	8x12 mm	No	No	1	Yes	Yes	Yes	2	Dual	No

For safe and proper use of products mentioned herein, please refer to the appropriate Instructions for Use.

Accessories Novaline
CE 0086 CE 0197

Baxter, AK 95 S, AK 96, AK 200 S, AK 200 ULTRA S, Artis,
Artis Physio, Integra and Novaline are trademarks of
Baxter International Inc. or its subsidiaries.

Baxter

서울시 중로구 중로 1
교보생명빌딩 10층 1001호
대표전화 : 02-6262-7100